Professional warewashing, cleaning and disinfection technology

Technical data sheet



TopClean M

Type: M2 Execution for: Oman

Technical data

Chemo-thermal disinfection 3-phase current: 3N PE 400V 50Hz Water connection: Cold water GiO MODULE (plinth, 120 mm)



Sample illustration

Rack capacity/h (theoretical)	10 / 6 / 5 racks/h
Programme cycle time *	360 / 540 / 720 s (3rd progamme running time with change of water
Rack dimension	500 x 500 mm to 540 x 500 mm
Entry height	420 mm
Dimensions (W x Hmin x D)	600 x 1340 x 600 mm
Electrical feeding cable	3-phase current 3N PE 400V 50Hz*
	nominal capacity: 12.4 kW
	nominal current: 19.6 A
Local fuse protection	25 A
Protection class of the machine	IP X5
Equipment	MIKE 2 control
	Infrared interface for wireless communication
	Leakage detector
	Soft start
	Boiler safety device
	Drain pump
	Automatic self-cleaning when tank is drained
	Back wall cladding
	Automatic door lock
	Compressed air connection for specially cleaned and de-oiled
	breathable air
	to EN 12021, 2-5 bar
	Integrated reverse osmosis (with connection and discharge set)
	connecting kit 1,6 m
	Integrated cupboard 400 mm
	Integrated cupboard, right hinged door
Fresh water line	For reverse osmosis: connecting set with stop valve,
	Pressure regulator, filter 10 µm with activated carbon
	Internal: air gap 'AA' in accordance with EN 1717 with booster pum
Fresh water supply	Minimum flow pressure 100 kPa / 1.0 bar in front of solenoid valve
	Maximum pressure 500 kPa / 5.0 bar in front of pressure reducer

Technical data sheet



Flow rate	total 4.3 l/min (at 15 °C inlet temperature and a flow pressure of 20 kPa / 2 bar)
	Permeate 2.0 I/min
	Concentrate 2.3 I/min
Fresh water thresholds for GiO MODULE	Temperature cold water min. 1 °C to max. 35 °C
	electrical conductibility < 1000 μS/cm
	Total hardness max. 28 / 5 mmol/I CaCO3
	Free of particles > 10 μm
	Iron < 0.10 mg/l
	Manganese < 0.04 mg/l
	Chlorine < 0.10 mg/l (free chlorine)
	Potassium permanganate < 10 mg/l
	Silicic acid < 10 mg/l
Final rinse water quantity	4.2 liters/cycle, variable
Boiler	Temperature: 50 °C
	Contents: 7.0 I
Wash tank	Filling: 11.0 I
	Temperature: 60 °C
Dosing of rinse aid	Hose pump (24 V) with time control
	and suction lance with empty level monitoring
Detergent dosage	Hose pump (24 V) with time control
	and suction lance with empty level monitoring
Material	Cladding: 1.4301
	Wash tank: 1.4301
	Boiler: 1.4404
Heat emission	for 6 programme cycles/h
	total: 0.6 kW
	perceptible: 0.4 kW
	latent: 0.2 kW
Ventilation flow rate	160 m³/h
Steam emission	0.3 kg/h
Emission sound pressure level at the workplace (LpA)	65 dB
Net / gross weight	102.5 kg / 118.5 kg (standard packaging)
Packaging dimensions (W x H x D)	700 x 1440 x 700 mm (standard packaging)

*Note:

Electrical equipment suitable for supply voltage: 3N PE 400 V 50 HZ (3N PE 380-415 V 50 Hz)